

PRIMEDIC™ Emergency cart



Special features:

- practical and modular design
- easy to move
- four swivel casters with two brakes
- prepared for PRIMEDIC™ wall fixation
- robust structure
- wall deflector
- one directional lock
- plastic cover white
- extremely good cost - performance ratio

General:

size	63 x 86 x 58 cm (W x H x D)
weight	55 kg
color	RAL 9002
Drawer Type 3D	21,8 x 14,5 x 47,5 cm (W x H x D)
Drawer Type 6D	52 x 14,5 x 47,5 cm (W x H x D)

Supply with Emergency Cart:

1	Emergency Cart	
3	drawers	
2	swivel casters	Ø 125 mm
2	swivel casters with brake	Ø 125 mm

Accessories available on demand. ZORYN MEDICAL CORP. WWW.ZORYN.COM E-MAIL: ZORYN@ATT.NET